

IMPORTANT SAFEGUARDS

Please read carefully before use in order to use the product in the correct manner and to prevent dangers or material damages in advance.

READ ALL INSTRUCTIONS BEFORE USING

DANGER:

1. Do not touch the power plug with wet hands under any circumstances.
2. Make sure that the product does not become wet.
3. Always unplug this product immediately after using.
4. Do not use while bathing.
5. Do not place or store product where it can fall or be pulled into a tub or sink.
6. Do not spill water or other liquid into the unit.
7. Do not reach for a product that has fallen into water.

WARNING:

1. The following persons should not use this product under any circumstances.
 - Person with acute diseases
 - Pregnant Woman
 - Person with malignant tumors
 - Person suffering from high fevers
 - Diabetes patient or person with infectious diseases
 - Person with skin allergies
 - Person with other unlisted diseases should consult with the doctor before use.
2. Do not use near volatile or flammable objects.
3. Do not use in humid areas.
4. Make sure the plastic wrapping material of the product is not accessible to children.
5. Try not to use the unit at temperatures of over 120°F if possible.
6. Do not use the product if you have cuts on your hand.
7. Do not put materials such as paper, fabric or vinyl in the hand entrance
(inside the unit).
8. A product should never be left unattended when plugged in.
9. Close supervision is necessary when this product is used by, on, or near children.
10. Use this product only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
11. Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the product to a service center for examination and repair.
12. Keep the cord away from heated surfaces.
13. Never block the air openings of the product or place it on a soft surface, such as a bed or couch, where the air openings may be blocked. Keep the air openings free of lint, hair, and the like.
14. Never use while sleeping or drowsy.
15. Never drop or insert any object into any opening or hose.

16. Do not use outdoors or operate where aerosol (spray) products are being used or where oxygen is being administered.
17. Connect this product to a properly grounded outlet only. See grounding instructions.
18. Unplug this product before filling. Fill(reservoir) with water only unless otherwise specified by manufacturer. Do not overfill (or specify filling instructions).

CAUTION:

1. Do not sleep with your hands are inserted into the product.
2. Do not touch the protective screen of the halogen lamp inside the main unit.
3. Try not to drop the product while moving it (damage may occur as the heating element is a loess magnetic board.)
4. Do not scratch the product with sharp objects.
5. Unplug the unit when it is not being used for long periods of time.
6. Do not plug in multiple products into an outlet at once.
7. If you feel any physical discomfort while using the unit, stop using immediately and consult with a doctor.
8. In cases of product malfunctions (unusual noises, smell, smoke, etc) unplug the power cord and contact your local center.
9. In case of use by a child, always use the unit under the supervision of an adult.
10. Do not allow infants to approach during use.
11. Do not block the air intake and take precautions so that objects (laundry, curtains, paper, etc) do not go in or block it.
12. Do not operate the unit with the air intake open, or open the intake during use.
13. Do not disassemble, repair, alter or modify the unit at your own discretion. Also, do not open the cover of the floor panel. (Not subject to After Service)
14. Do not clean the floor loess board with water or detergents.
15. Do not wash main unit with water or detergents.
16. Be careful in moving the unit with your hand inside it as the unit may drop or electrical hazards may occur.
17. Be careful not to look directly at the infrared emission source.
18. Do not use for longer than 60 minutes during a single session.
19. If you apply substances such as lotion to your hand before use the loess floor board will get stained easily.
20. Refrain from placing your hand past the floor board boundary line. (danger of burn)
21. Do not cook on the floor board.

SAVE THESE INSTRUCTIONS

1. Functions of the product

1) Functions of Skinfive

①Loess forms a multi-layered and multi-punctured structure with large surface areas that allow it to absorb and save large amounts of far-infrared rays within to release these large amounts of far-infrared rays when heated. Skinfive uses this characteristic of loess to heat the loess floor board and loess emission unit in order to emit large amounts of infrared rays. (Korea Institute of Construction Materials test grade NO- FIR482, 483)

Also, we have proven that our product emits large amounts of far-infrared rays by acquiring the 'Won-mark', the symbol for far-infrared rays, from the Institute of Construction Materials.

② Loess emits large quantities of natural 'negative-ions', also known as vitamins of the air.

(Korea Institute of Construction Materials test grade NO- FIR349, 350)

③ Loess has a natural smell removing quality. Results of tests carried out in respect to the actual smell removing capabilities of Skinfive's loess floor board and loess emission unit show that our product is effective in removing the smell ammonia, the main cause of bad smells, and also formaldehyde, the characteristic toxin of Sick House Syndrome, thus creating a healthy atmosphere.

(Korea Institute of Construction Materials test grade NO- FIR391, 394)

2) Skinfive Thermo Treatment(heat treatment) Function

① The loess floor board, baked at 2246°F has a thermo treatment(heat treatment) function that helps relieve muscle tensions.

② The Skinfive floor board possesses the natural curves and the pressure protrusions offer heat treatment and pressure massage functions.

3)Skinfive Infrared Functions

① The light and heat of the infrared halogen lamp radiate on the loess emission unit, irradiates infrared rays and heat to penetrate deep into the skin to relieve muscle tensions and pains.

4) How to take care of and manage the unit after use

① Remove the power plug from the outlet when putting the unit aside after use.

② When not using the unit for long periods of time, take precautions so that the unit is not damaged due to dust, moisture, mildew etc.

③ Preserve in a dry area at room temperatures and humidity, away from direct sunlight.

Cleaning the Unit

- When cleaning the outside of the unit, unplug the power plug and wipe lightly using a soft cloth.

Clean the filter located at the back of the fan cover on a regular basis so as to provide the unit with clean air.

Wipe the loess floor board lightly with a soft, dry cloth.

Using the Skinfive effectively

1) Recommended sessions of use of Skinfive per day (Temperature, Time)

Generally, it is recommended that you use 1~2 times a day, in the morning and evening, for 20 minutes at a time at temperatures between 104°F~111°F in order to experience positive effects on your skin.

Users who wish to experience relief of muscle and hand pains should use 2~3 times a day for 20~40 minutes at a time at temperatures of 113°F~118°F

2) Correct Usage

- The appropriate use of Skinfive differs according to the condition of the hand of the user, seasonal factors, and the purpose of use. First of all, wash your hands thoroughly before using this product and stick to the recommended sessions and times of use per day.
- Refrain from using the unit at temperatures above 113°F. If you feel that your hands are too hot at any point during usage, remove your hands and use again a little while later.
- When using the unit for long periods of time, remove your hands from the unit from time to time and re-insert hands a while later.
- Skinfive will be even more effective if you apply moisturizer to your hands after use.

3) Where to use Skinfive

Skinfive is easy to carry around and thus may be easily used in all areas including on your desk, table, or dresser. You can place the unit on your lap and insert your hands to operate it and enjoy watching TV, or carrying on conversations, and you can also operate the unit on one hand while using the other to take care of other tasks (making phone calls, writing etc.).

Self-Diagnosis

1. Steps to take before requesting after service

Categories	Steps to take
When the unit does not switch on	<ul style="list-style-type: none">① Check if the power cord is plugged into an outlet② Plug in other products to the outlet to see if electricity is supplied
When the brightness of the halogen lamp is unstable	<ul style="list-style-type: none">① This is normal as the brightness of the lamp is controlled automatically when the unit reaches the temperature set by the user.② When the halogen lamp won't turn on press the power button again, and contact a nearby A/S center if the problem persists
When wind is not properly induced into the unit	<ul style="list-style-type: none">① Check the air intake located at the back of the unit for objects that may interfere with air intake.② Turn the air intake at the back of the unit to open, clean the filter softly with water, dry, and insert the filter back.
When you smell something unusual	<ul style="list-style-type: none">① The unit may smell a little when first purchased as it is a new product.② When you smell dust being sucked in clean the filter which is located at the back of the product with water, dry, and put it back in.③ If you constantly smell a something burning unplug the power cord and contact a nearby A/S center.

<p>When the unit does not operate even after you press the ON button and insert your hands into the unit</p>	<p>① Check to see if your hands are placed accurately over the sensors located at the bottom of the hand entrances. ② If that is not the problem turn the unit OFF and then back ON, insert your hands into the hand entrances and place the directly over the sensors. ③ If the problem persists contact a nearby A/S center.</p>
--	--